

ENDOSCOPE TRANSPORT BAG

The CatchAll Transport Pad is a convenient way to transport a soiled endoscope from the procedure to reprocessing room. When securely closed, the pad cushions and protects the scope while providing a sealed transport container. The pad is simple to set up; includes are sizeable absorbent material pad and thoroughly covers the cart surface. After use, cinch up the drawstrings and the pad converts to a sealed, cushioned container. The thick absorbent material helps both absorb any fluid and protect the scope during transportation.





Specifications:

Transport Pad		
REF	Description	Package Units
CTP-01025	CatchAll™Transport Pad	100/Box