



CK SURGITECH

HELICOBACTER PYLORI BREATH TESTING

What is H. Pylori?

A spiral-shaped, gram-negative, curved rod bacterium. 50% of the world's population are infected.

Who Discovered H. Pylori?





Robin Warren and
Barry Marshall




How are people infected
with H. Pylori?

Oral-Oral, Fecal-Oral
or Gastro-Oral



Stomach Ulcer Formation



CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214
P: 07 5574 5573 E: sales@cksurgitech.com



CK SURGITECH

HELICOBACTER PYLORI BREATH TESTING

Sensitivity 95%-99%, specificity 95%-99%.

Results obtained on site, within 4 minutes of sampling

Easy to use, near patient testing, no need to send samples for analysis.

Non-invasive procedure. Swallow Urea capsule and take the breath sample

No cross infection; Detection process doesn't contain any hazardous substances.

No need for complex and expensive test equipment, analysis made on-site.



Take a Urea Capsule with water



Sit for 15 minutes



Blow into the collection card for 3-5 minutes



Test the collection card with the analyzer to obtain the result



CK Surgitech Pty Ltd

6/14 Technology Drive, Arundel
Queensland Australia 4214
P: 07 5574 5573 E: sales@cksurgitech.com